





C-LOCK



INSTRUCTIONS FOR USE (USE ASEPTIC TECHNIQUE)

- 01** Prior to initiation of a venous access session, aspirate and discard the solution that was instilled in the dialysis catheter during the previous session.
- 02** At the completion of each venous access session, flush the lumens of the catheter with 5mL of sterile Tribasic Sodium Chloride solution to remove any blood, medication or any other substance remaining in the catheter.
- 03** Aseptically open package containing C-LOCK prefilled syringe(s). Inspect the contents of the syringe(s) according to the provided instructions for use.
- 04** Hold the syringe with the cap end up and carefully twist the cap off.

- 05** To remove any air bubbles, point luer end of syringe up and gently tap the syringe with fingers until air bubbles rise up to top of the syringe.
- 06** Keep luer end of syringe pointing up and slowly push the plunger up to force air bubbles out of the syringe.
- 07** Properly identify priming volume (A) and with the luer end of syringe pointing up, slowly push plunger to the line on the syringe that matches the priming volume (B) (prescribed by the catheter manufacturer).


- 08** Connect the syringe to the catheter. Use only with compatible luer lock connectors.

- 09** Slowly, over a 5 to 10 second period, instill the exact priming volume, as specified by the catheter manufacturer, of C-LOCK solution into each catheter lumen, as a locking agent.
- 10** Remove and discard unused portions and empty syringes according to your biohazard waste policy.



D.B.M. S.r.l. a socio unico
Via dell'Industria, 2/8
23030 Tovo di Sant'Agata (SO)
ITALY

D.B.M. SALINEFLUSH is manufactured in Italy in accordance with the ISO 13485 Standard to guarantee quality and reliability.



FS521104

ARTG #497145
ARTG #497185

TekMed™
1300 720 727
www.tekmed.com.au

TekMed NZ™
0800 456 217
www.tekmed.co.nz

