

SPLIT STREAM®

FEATURES & BENEFITS

- Retrograde split tip catheter allows for more accurate placement.
- Available in 14F and 16F in multiple lengths for all patient needs.
- Removable hub aids in patient comfort.
- Over-the-wire version permits a “sheathless” technique, especially useful for left-sided procedures.

Catalog #		
SST24E.	14F x 24CM RETROGRADE SPLIT STREAM®, BASIC SET	5/Box
SST28E.	14F x 28CM RETROGRADE SPLIT STREAM®, BASIC SET	5/Box
SST32E.	14F x 32CM RETROGRADE SPLIT STREAM®, BASIC SET	5/Box
SST36E.	14F x 36CM RETROGRADE SPLIT STREAM®, BASIC SET	5/Box
SST40E.	14F x 40CM RETROGRADE SPLIT STREAM®, BASIC SET	5/Box
SST2416E.	16F x 24CM RETROGRADE SPLIT STREAM®, BASIC SET	5/Box
SST2816E.	16F x 28CM RETROGRADE SPLIT STREAM®, BASIC SET	5/Box
SST3216E.	16F x 32CM RETROGRADE SPLIT STREAM®, BASIC SET	5/Box
SST3616E.	16F x 36CM RETROGRADE SPLIT STREAM®, BASIC SET	5/Box
SST4016E.	16F x 40CM RETROGRADE SPLIT STREAM®, BASIC SET	5/Box

Catalog #	Split Stream® without Sideholes	
SST24WOE.	14F x 24CM SPLIT STREAM® W/O SIDEHOLES, BASIC SET	5/Box
SST28WOE.	14F x 28CM SPLIT STREAM® W/O SIDEHOLES, BASIC SET	5/Box
SST32WOE.	14F x 32CM SPLIT STREAM® W/O SIDEHOLES, BASIC SET	5/Box
SST36WOE.	14F x 36CM SPLIT STREAM® W/O SIDEHOLES, BASIC SET	5/Box
SST40WOE.	14F x 40CM SPLIT STREAM® W/O SIDEHOLES, BASIC SET	5/Box


SET CONTAINS:

- (1) Catheter (1) Guidewire (2) Dilators (1) Arterial Extension Set (Red) (1) Venous Extension Set (Blue) (1) Valved Peelable Introducer (4) Clamps (1) Removable Hub (2) Tunnelers (1) Lumen Clamp (2) Injection Caps (1) Introducer Needle

Catalog #	Split Stream® OTW	
SST24SE.	14F x 24CM SPLIT STREAM®, BASIC SET	5/Box
SST28SE.	14F x 28CM SPLIT STREAM®, BASIC SET	5/Box
SST32SE.	14F x 32CM SPLIT STREAM®, BASIC SET	5/Box
SST36SE.	14F x 36CM SPLIT STREAM®, BASIC SET	5/Box

OTW SET CONTAINS:

- (1) Catheter w/ Stylet (1) 120cm Teflon Coated Guidewire (1) 70cm J-Flex Guidewire (2) Dilators (1) Arterial Extension Set (Red) (1) Venous Extension Set (Blue) (1) Valved Peelable Introducer (4) Clamps (1) Removable Hub (2) Tunnelers (1) Lumen Clamp (2) Injection Caps (1) Introducer Needle



TekMed

1300 720 727

12 Bellevue Crescent
Preston Victoria 3072 Australia

T (03) 9487 9999 F (03) 9484 8792
info@tekmed.com.au • www.tekmed.com.au