

## FREE FLOW™ ORDERING INFORMATION



**F  
R  
E  
E  
F  
L  
O  
W™**

### Long-Term Free Flow

15.5FR Straight	Sets (5/Box)
15.5FR x 24cm Straight Long-Term Free Flow Catheter with Stylet Set	MFFL1524S
15.5FR x 28cm Straight Long-Term Free Flow Catheter with Stylet Set	MFFL1528S
15.5FR x 32cm Straight Long-Term Free Flow Catheter with Stylet Set	MFFL1532S
15.5FR x 36cm Straight Long-Term Free Flow Catheter with Stylet Set	MFFL1536S

15.5FR Pre-Curved	Sets (5/Box)
15.5FR x 24cm Pre-Curved Long-Term Free Flow Catheter with Stylet Set	MFFL1524PC
15.5FR x 28cm Pre-Curved Long-Term Free Flow Catheter with Stylet Set	MFFL1528PC
15.5FR x 32cm Pre-Curved Long-Term Free Flow Catheter with Stylet Set	MFFL1532PC
15.5FR x 36cm Pre-Curved Long-Term Free Flow Catheter with Stylet Set	MFFL1536PC

Set contains: (1) Catheter with Stylet, (1) Grasper Tunneler, (1) Peel-Away Sheath, (2) Vessel Dilators, (1) J-Flex Guidewire, (1) Introducer Needle, (2) Injection Caps

### Short-Term Free Flow

12.5FR Straight	Sets (5/Box)
12.5FR X 15cm Straight Short-Term Free Flow with Stylet Catheter Set	MFFS1215S
12.5FR X 20cm Straight Short-Term Free Flow with Stylet Catheter Set	MFFS1220S
12.5FR X 24cm Straight Short-Term Free Flow with Stylet Catheter Set	MFFS1224S

15.5FR Straight	Sets (5/Box)
15.5FR X 15cm Straight Short-Term Free Flow with Stylet Catheter Set	MFFS1515S
15.5FR X 20cm Straight Short-Term Free Flow with Stylet Catheter Set	MFFS1520S
15.5FR X 24cm Straight Short-Term Free Flow with Stylet Catheter Set	MFFS1524S

12.5FR Raulerson IJ	Sets (5/Box)
12.5FR X 12cm IJ Short-Term Free Flow with Stylet Catheter Set	MFFS1212IJ
12.5FR X 15cm IJ Short-Term Free Flow with Stylet Catheter Set	MFFS1215IJ
12.5FR X 20cm IJ Short-Term Free Flow with Stylet Catheter Set	MFFS1220IJ

15.5FR Raulerson IJ	Sets (5/Box)
15.5FR X 12cm IJ Short-Term Free Flow with Stylet Catheter Set	MFFS1512IJ
15.5FR X 15cm IJ Short-Term Free Flow with Stylet Catheter Set	MFFS1515IJ
15.5FR X 20cm IJ Short-Term Free Flow with Stylet Catheter Set	MFFS1520IJ

12.5FR IJ with 2 Wings	Sets (5/Box)
12.5FR x 12cm IJ Short-Term Free Flow with Stylet and 2 Suture Wings	MFFS1212IJ-2
12.5FR x 15cm IJ Short-Term Free Flow with Stylet and 2 Suture Wings	MFFS1215IJ-2
12.5FR x 20cm IJ Short-Term Free Flow with Stylet and 2 Suture Wings	MFFS1220IJ-2

15.5FR IJ with 2 Wings	Sets (5/Box)
15.5FR x 12cm IJ Short-Term Free Flow with Stylet and 2 Suture Wings	MFFS1512IJ-2
15.5FR x 15cm IJ Short-Term Free Flow with Stylet and 2 Suture Wings	MFFS1515IJ-2
15.5FR x 20cm IJ Short-Term Free Flow with Stylet and 2 Suture Wings	MFFS1520IJ-2

Set contains: (1) Catheter with Stylet, (3) Vessel Dilators, (1) J-Flex Guidewire, (1) Introducer Needle, (2) Injection Caps



**No Sideholes / Open Tips**

 **medCOMP**  
 1499 Delp Drive  
 Harleysville, PA 19438  
 Tel: 215-256-4201  
 Fax: 215-256-1787  
 www.medcompnet.com

**TekMed**  
**T (03) 9487 9999**  
 12 Bellevue Crescent  
 Preston Victoria 3072 Australia  
**F (03) 9484 8792**  
 info@tekmed.com.au • www.tekmed.com.au

PN 2476 REV. B 9/08

www.medcompnet.com



# FREE-FLOW™

Short & Long Term Hemodialysis Catheter

*As innovation continues superior products come to life.*

Medcomp® presents its premier *open tip* dialysis catheter, the Free Flow™. Through the use of leading technology and unprecedented engineering the industry's first true 360° coaxial flow design has come to life, eliminating the need for side holes and reversing of catheter lines due to arterial insufficiency. The *lack of sideholes* preserves catheter locking solution concentration and minimizes locking solution leaching. This unique design increases the overall patency and performance of the catheter.

## Rotating Suture Wing

Rotating suture wing provides external anchoring.



## Crack Resistant Luers

The ergonomical design ensures greater resistance to cracking or separation.

## Pre-Loaded Stylet

Allows for standard or modified Seldinger introduction. Once removed, renders the venous lumen totally open, unlike tapered tipped catheters.

## True 360° Flow

- Proven Coaxial Design
- Reduces Arterial Insufficiencies



## No Sideholes

Reduces loss of DuraLock-C™ and other catheter lock solutions.

## Open Tips

Unobstructed uptake and return reduce pressure and maximize flow.

## Raulerson IJ



## Soft Tip Stylet (Short Term)

Integral soft tip reduces the risk of vessel damage.



FREE-FLOW™

