



At 15.5F, the T-3 CT catheter is the largest triple lumen power injectable catheter available. Its large lumen size enables greater flow rates, which improve the administration of extracorporeal therapies including dialysis, hemofiltration, apheresis and plasmapheresis. A power injectable central lumen expands the device's capability by enabling infusion, power



injection of contrast media and central venous pressure monitoring. Manufactured from radiopaque thermosensitive material, the T-3 CT provides increased patient comfort, excellent resistance to kinking and improved biocompatibility. The T-3 CT marks another step forward in Medcomp's continued goal of improving clinical access and patient's quality of life.



## Double Lumen Reliability, CT Injectable Triple Lumen Versatility

- 1 Triple Lumen**  
Two dedicated dialysis lumens coupled with a power injectable central lumen provides convenient access for not only dialysis but infusion, power injection of contrast media and withdrawal needs, enabling multiple therapies through one access site.
- 2 Radiopaque Polyurethane**  
The radiopaque polyurethane catheter material provides rigidity for insertion, kink resistance to ensure adequate treatment, and thermosensitivity to provide increased patient comfort.
- 3 Available in 15, 20, 24, 28 & 32cm Insertion Length**  
Ranging sizes from 15 to 32cm, the T-3 CT catheter accommodates various insertion sites and individual patient anatomy to achieve the most adequate treatments possible.
- 4 Step Tip Design**  
With an open Step Tip Design and a 90-degree arterial lumen, the T-3 CT catheter provides optimal flow rates, reduced arterial insufficiency and improves the ability to hold locking solutions, maintaining catheter patency.
- 5 15.5F Lumen**  
The T-3 CT's large lumen size enables greater flow rates, which improve the administration of extracorporeal therapies including dialysis, hemofiltration, apheresis and plasmapheresis.
- 6 Rotatable Suture Wing**  
Rotatable Suture Wing provides increased patient comfort and device securement.
- 7 Isoplast Luers**  
The isoplast luers have an outstanding practical toughness, dimensional stability and chemical resistance to provide you and your patients with a safe and reliable dialysis connection for the duration of treatment.
- 8 Priming Volumes**  
Clearly labeled priming volumes ensure a proper and adequate instillation of locking solution to maintain optimal catheter patency.

[medcompnet.com/T-3](http://medcompnet.com/T-3)



# THE LARGEST TRIPLE LUMEN POWER INJECTABLE CATHETER AVAILABLE.



	CATHETER ONLY	UoM
MC013541	15.5F X 15CM T-3 CT, STRAIGHT EXTENSION	5/BOX
MC013542	15.5F X 20CM T-3 CT, STRAIGHT EXTENSION	5/BOX
MC013543	15.5F X 24CM T-3 CT, STRAIGHT EXTENSION	5/BOX
MC013544	15.5F X 28CM T-3 CT, STRAIGHT EXTENSION	5/BOX
MC013545	15.5F X 32CM T-3 CT, STRAIGHT EXTENSION	5/BOX

#### CATHETER ONLY CONTAINS:

- (1) Catheter
- (3) Dilators
- (1) Valved Introducer
- (1) Patient ID Card
- (3) Injection Caps

	BASIC SET	UoM
MC013501	15.5F X 15CM T-3 CT, STRAIGHT EXTENSION	5/BOX
MC013502	15.5F X 20CM T-3 CT, STRAIGHT EXTENSION	5/BOX
MC013503	15.5F X 24CM T-3 CT, STRAIGHT EXTENSION	5/BOX
MC013504	15.5F X 28CM T-3 CT, STRAIGHT EXTENSION	5/BOX
MC013505	15.5F X 32CM T-3 CT, STRAIGHT EXTENSION	5/BOX

#### BASIC SET CONTAINS:

- (1) Catheter
- (3) Dilators
- (1) Scalpel
- (1) Valved Introducer
- (1) Guidewire
- (1) Adhesive Wound Dressing
- (1) 2-0 Silk Suture
- (1) Patient ID Card
- (3) Injection Caps
- (1) Introducer Needle


**medCOMP**  
 1499 Delp Drive  
 Harleysville, PA 19438  
 Tel: 215-256-4201  
 Fax: 215-256-1787  
[www.medcompnet.com](http://www.medcompnet.com)

[medcompnet.com/T-3](http://medcompnet.com/T-3)



**TekMed™**

**1300 720 727**

12 Bellevue Crescent Preston Victoria 3072 Australia

T (03) 9487 9999 F (03) 9484 8792

info@tekmed.com.au www.tekmed.com.au



FS521104

PN2563 Rev. A 11/13