

PERIPHERALLY INSERTED
CENTRAL CATHETER (PICC)
VASCU-PICC® 2.6F



One size does not fit all.

1.9F SINGLE LUMEN PICC

Catalog Code	Description	UOM	QTY
MR17011101	1.9F X 50CM SINGLE LUMEN PICC CATHETER W/ STYLET	EA	5/Box
MR17011121	1.9F X 20CM SINGLE LUMEN PICC CATHETER W/ STYLET	EA	5/Box

Tray Components:
 Vascu Picc w/ Preloaded Stylet (1) 10cc Syringe (1) 2F (22ga) Tearaway
 (1) Latex Free Tourniquet (1) Needle Free Valve (1) Suture Wing
 (10) 2 x 2 Gauze (1) 24" Tape Measure (1) Patient Information Booklet
 (1) Catheter Securement Device

Catalog Code	Description	UOM	QTY
VP1.9S20-NS	1.9F X 20CM SINGLE LUMEN PICC CATHETER W/O STYLET	EA	5/Box
VP1.9S50-NS	1.9F X 50CM SINGLE LUMEN PICC CATHETER W/O STYLET	EA	5/Box

Tray Components:
 Vascu Picc Catheter (1) 10cc Syringe (1) 2F (22ga) Tearaway
 (1) Latex Free Tourniquet (1) Needle Free Valve (1) Suture Wing
 (10) 2 x 2 Gauze (1) 24" Tape Measure (1) Patient Information Booklet
 (1) Catheter Securement Device

2.6F DUAL LUMEN PICC

Catalog Code	Description	UOM	Qty
MR17012601	2.6F X 50CM DUAL LUMEN PICC W/ STYLET	EA	5/Box
MR17012621	2.6F X 20CM DUAL LUMEN PICC W/ STYLET	EA	5/Box

Tray Components:
 Vascu Picc w/ Preloaded Stylet (1) 10cc Syringe (1) Catheter Securement Device
 (1) Latex Free Tourniquet (2) Needle Free Valves (1) Patient Information Booklet
 (10) 2 x 2 Gauze (1) 3F (19ga) Tearaway (1) Patient Chart Sticker
 (1) 24" Tape Measure

Catalog Code	Description	UOM	QTY
MR17012628	2.6F X 20CM DUAL LUMEN PICC CATHETER W/ CUFF AND STYLET	EA	5/Box
MR17012608	2.6F X 50CM DUAL LUMEN PICC CATHETER W/ CUFF AND STYLET	EA	5/Box

Tray Components:
 Vascu Picc w/ Preloaded Stylet (2) 10cc Syringe (1) Catheter Securement Device
 (1) Latex Free Tourniquet (2) Needle Free Valves (1) Patient Information Booklet
 (10) 2 x 2 Gauze (1) 3F (19ga) Tearaway (2) Patient Chart Sticker
 (1) 24" Tape Measure

PERIPHERALLY INSERTED
CENTRAL CATHETER (PICC)
VASCU-PICC® 2.6F



TekMed
1300 720 727
 12 Bellevue Crescent
 Preston Victoria 3072 Australia
 T (03) 9487 9999 F (03) 9484 8792
 info@tekmed.com.au • www.tekmed.com.au

medcompnet.com/2.6F

DUAL LUMEN VASCU-PICC® 2.6F

SMALL but POWERFUL

Approved for infusion therapy

From the smallest to the largest — and now everyone in-between — Medcomp® has a catheter size for any pediatric patient. VASCU-PICC® 2.6F

is the latest addition to our PICC product line and features a dual lumen configuration, soft polyurethane body, and high flow rates in a small French size. Since smaller veins require greater precision, it only makes sense to offer you more choices in PICC sizes.

Comfort

Comfort is a special challenge with your pediatric patients. The smaller French size features a dual lumen which means fewer sticks, fewer tears.

Versatility

Since many pediatric patients require a combination of nutrition, fluids and/or medicine, the dual lumen configuration enables one catheter to serve multiple functions.

Perfect Fit

Medcomp® offers a full-line of pediatric PICC sizes to fit any pediatric patient. The new 2.6F fits nicely between existing 1.9F and 3F catheters, and makes catheter insertion easier.

Multiple Function Luers

Eliminates separate access lines and decreases vessel trauma

Dual Lumen Design

- Combined access to one site reduces chance of infection, providing better patient outcomes and increased parent satisfaction
- Offers the ability to infuse incompatible drugs




Especially with your smallest patients.

Thermosensitive Material
Softens at body temperature

Cuff
Promotes tissue in-growth, securing the catheter and reducing infection

**Easy-to-read
Depth Marks**



Radiopaque
Tip identifiable under X-Ray

Polyurethane Material

- Provides strength, versatility, and high flow rates in a smaller French size
- Higher flow rates and less chance of occlusion
- Allows greater control during insertion

2.6F Catheter Size 20cm & 50cm Lengths

Smaller catheter while maintaining dual access and reduced trauma in smaller vessels

