

FREE FLOW™ ORDERING INFORMATION



**F
R
E
E
F
L
O
W™**

Long-Term Free Flow

15.5FR Straight		Sets (5/Box)
15.5FR x 24cm Straight Long-Term Free Flow Catheter with Stylet Set	MFFL1524S	
15.5FR x 28cm Straight Long-Term Free Flow Catheter with Stylet Set	MFFL1528S	
15.5FR x 32cm Straight Long-Term Free Flow Catheter with Stylet Set	MFFL1532S	
15.5FR x 36cm Straight Long-Term Free Flow Catheter with Stylet Set	MFFL1536S	

15.5FR Pre-Curved		Sets (5/Box)
15.5FR x 24cm Pre-Curved Long-Term Free Flow Catheter with Stylet Set	MFFL1524PC	
15.5FR x 28cm Pre-Curved Long-Term Free Flow Catheter with Stylet Set	MFFL1528PC	
15.5FR x 32cm Pre-Curved Long-Term Free Flow Catheter with Stylet Set	MFFL1532PC	
15.5FR x 36cm Pre-Curved Long-Term Free Flow Catheter with Stylet Set	MFFL1536PC	

Set contains: (1) Catheter with Stylet, (1) Grasper Tunneler, (1) Peel-Away Sheath, (2) Vessel Dilators, (1) J-Flex Guidewire, (1) Introducer Needle, (2) Injection Caps

Short-Term Free Flow

12.5FR Straight		Sets (5/Box)
12.5FR X 15cm Straight Short-Term Free Flow with Stylet Catheter Set	MFFS1215S	
12.5FR X 20cm Straight Short-Term Free Flow with Stylet Catheter Set	MFFS1220S	
12.5FR X 24cm Straight Short-Term Free Flow with Stylet Catheter Set	MFFS1224S	

15.5FR Straight		Sets (5/Box)
15.5FR X 15cm Straight Short-Term Free Flow with Stylet Catheter Set	MFFS1515S	
15.5FR X 20cm Straight Short-Term Free Flow with Stylet Catheter Set	MFFS1520S	
15.5FR X 24cm Straight Short-Term Free Flow with Stylet Catheter Set	MFFS1524S	

12.5FR Raulerson IJ		Sets (5/Box)
12.5FR X 12cm IJ Short-Term Free Flow with Stylet Catheter Set	MFFS1212IJ	
12.5FR X 15cm IJ Short-Term Free Flow with Stylet Catheter Set	MFFS1215IJ	
12.5FR X 20cm IJ Short-Term Free Flow with Stylet Catheter Set	MFFS1220IJ	

15.5FR Raulerson IJ		Sets (5/Box)
15.5FR X 12cm IJ Short-Term Free Flow with Stylet Catheter Set	MFFS1512IJ	
15.5FR X 15cm IJ Short-Term Free Flow with Stylet Catheter Set	MFFS1515IJ	
15.5FR X 20cm IJ Short-Term Free Flow with Stylet Catheter Set	MFFS1520IJ	

12.5FR IJ with 2 Wings		Sets (5/Box)
12.5FR x 12cm IJ Short-Term Free Flow with Stylet and 2 Suture Wings	MFFS1212IJ-2	
12.5FR x 15cm IJ Short-Term Free Flow with Stylet and 2 Suture Wings	MFFS1215IJ-2	
12.5FR x 20cm IJ Short-Term Free Flow with Stylet and 2 Suture Wings	MFFS1220IJ-2	

15.5FR IJ with 2 Wings		Sets (5/Box)
15.5FR x 12cm IJ Short-Term Free Flow with Stylet and 2 Suture Wings	MFFS1512IJ-2	
15.5FR x 15cm IJ Short-Term Free Flow with Stylet and 2 Suture Wings	MFFS1515IJ-2	
15.5FR x 20cm IJ Short-Term Free Flow with Stylet and 2 Suture Wings	MFFS1520IJ-2	

Set contains: (1) Catheter with Stylet, (3) Vessel Dilators, (1) J-Flex Guidewire, (1) Introducer Needle, (2) Injection Caps



No Sideholes / Open Tips

medCOMP
 1499 Delp Drive
 Harleysville, PA 19438
 Tel: 215-256-4201
 Fax: 215-256-1787
 www.medcompnet.com

TekMed
 1300 720 727
 12 Bellevue Crescent Preston Victoria 3072 Australia
 T (03) 9487 9999 F (03) 9484 8792
 info@tekmed.com.au www.tekmed.com.au

PN 2476 REV. B 9/08

www.medcompnet.com



FREE-FLOW™

Short & Long Term Hemodialysis Catheter

As innovation continues superior products come to life.

Medcomp® presents its premier *open tip* dialysis catheter, the Free Flow™. Through the use of leading technology and unprecedented engineering the industry's first true 360° coaxial flow design has come to life, eliminating the need for side holes and reversing of catheter lines due to arterial insufficiency. The *lack of sideholes* preserves catheter locking solution concentration and minimizes locking solution leaching. This unique design increases the overall patency and performance of the catheter.

Rotating Suture Wing

Rotating suture wing provides external anchoring.



Crack Resistant Luers

The ergonomical design ensures greater resistance to cracking or separation.



Raulerson IJ

Pre-Loaded Stylet

Allows for standard or modified Seldinger introduction. Once removed, renders the venous lumen totally open, unlike tapered tipped catheters.

True 360° Flow

- Proven Coaxial Design
- Reduces Arterial Insufficiencies



No Sideholes

Reduces loss of DuraLock-C™ and other catheter lock solutions.

Open Tips

Unobstructed uptake and return reduce pressure and maximize flow.

Soft Tip Stylet (Short Term)

Integral soft tip reduces the risk of vessel damage.



FREE-FLOW™

