

ORDERING INFORMATION

Hemo-Cath® LT

| CATALOG # | BASIC SET | UOM |
|-----------|----------------------------------|-------|
| SL18P | 8F X 18CM SILICONE HEMO-CATH® LT | 5/BOX |
| SL24P | 8F X 24CM SILICONE HEMO-CATH® LT | 5/BOX |

Set Includes: (1) Catheter (1) Introducer Needle (1) Adhesive Dressing (1) J/Flex Guidewire (1) Tunneler (1) Hemo-Cath® Hub Clip (2) Injection Caps (1) 10cc Syringe (1) Peelable Introducer

| CATALOG # | BASIC SET | UOM |
|-----------|-----------------------------------|-------|
| MC101041P | 10F X 15CM SILICONE HEMO-CATH® LT | 5/BOX |
| MC101042P | 10F X 18CM SILICONE HEMO-CATH® LT | 5/BOX |
| MC101043P | 10F X 20CM SILICONE HEMO-CATH® LT | 5/BOX |
| MC101044P | 10F X 22CM SILICONE HEMO-CATH® LT | 5/BOX |
| MC101045P | 10F X 24CM SILICONE HEMO-CATH® LT | 5/BOX |

Set Includes: (1) Catheter (1) Introducer Needle (1) Adhesive Dressing (1) J/Flex Guidewire (1) Tunneler (1) Hemo-Cath® Hub Clip (2) Injection Caps (1) 10cc Syringe (1) Peelable Introducer

| CATALOG # | BASIC SET | UOM |
|-----------|-------------------------------------|-------|
| MC101241 | 12.5F X 15CM SILICONE HEMO-CATH® LT | 5/BOX |
| MC101242 | 12.5F X 18CM SILICONE HEMO-CATH® LT | 5/BOX |
| MC101243 | 12.5F X 24CM SILICONE HEMO-CATH® LT | 5/BOX |

Set Includes: (1) Catheter (1) Introducer Needle (1) Adhesive Dressing (1) J/Flex Guidewire (1) Tunneler (1) Hemo-Cath® Hub Clip (1) Scalpel (1) Peelable Introducer (2) Injection Caps

| CATALOG # | STRAIGHT CATHETER BASIC SET | UOM |
|-----------|-------------------------------------|-------|
| SL28E. | 12.5F x 28CM SILICONE HEMO-CATH® LT | 5/BOX |
| SL32E. | 12.5F x 32CM SILICONE HEMO-CATH® LT | 5/BOX |

| CATALOG # | PRE-CURVED CATHETER BASIC SET | UOM |
|-----------|--|-------|
| SL28PCE. | 12.5F X 28CM PC SILICONE HEMO-CATH® LT | 5/BOX |
| SL32PCE. | 12.5F X 32CM PC SILICONE HEMO-CATH® LT | 5/BOX |

Set Includes: (1) Catheter (2) Injection Caps (1) Peelable Introducer (1) Introducer Needle (1) J/Flex Guidewire (1) Tunneler (1) Hemo-Cath® Hub Clip (*Straight Catheter Only*)

Hemo-Cath® ST

| CATALOG # | BASIC SET | UOM |
|-----------|----------------------------------|-------|
| SL12P | 8F X 12CM SILICONE HEMO-CATH® ST | 5/BOX |

Set Includes: (1) Catheter (2) Injection Caps (1) J/Flex Guidewire (1) Vessel Dilator (1) Adhesive Dressing (1) Hemo-Cath® Hub Clip (1) 10cc Syringe (1) Introducer Needle

| CATALOG # | BASIC SET | UOM |
|-----------|-------------------------------------|--------|
| SL15E | 11.5F X 15CM SILICONE HEMO-CATH® ST | 10/BOX |
| SL20E | 11.5F X 20CM SILICONE HEMO-CATH® ST | 10/BOX |
| SL24E | 11.5F X 24CM SILICONE HEMO-CATH® ST | 10/BOX |
| SDL136E | 13.5F X 15CM SILICONE HEMO-CATH® ST | 10/BOX |
| SDL138E | 13.5F X 20CM SILICONE HEMO-CATH® ST | 10/BOX |
| SDL139E | 13.5F X 24CM SILICONE HEMO-CATH® ST | 10/BOX |
| SDL1311E | 13.5F X 28CM SILICONE HEMO-CATH® ST | 5/BOX |
| SDL1314E | 13.5F X 35CM SILICONE HEMO-CATH® ST | 5/BOX |

Set Includes: (1) Catheter with Stylet (1) Vessel Dilator (1) Introducer Needle (2) Injection Caps (1) J/Flex Guidewire (1) Hemo-Cath® Hub Clip

HEMO-CATH®
(LT & ST)



medcompnet.com

medCOMP

TekMed
1300 720 727
www.tekmed.com.au

TekMed NZ
0800 456 217
www.tekmed.co.nz



FS521104



medCOMP

1499 Delp Dr.
Harleysville, PA 19438
P: 215.256.4201 F: 215.256.1787
www.medcompnet.com

PN2541 Rev. C 4/17



TekMed™
1300 720 727
www.tekmed.com.au

TekMed NZ™
0800 456 217
www.tekmed.co.nz



Our Hemo-Cath® ST provides the ideal catheter, in comparison to polyurethane, for short term access:

| PUBLISHED RESULTS | | |
|--|-----------------------|-------------------|
| Temporary Access: Silicone vs. Polyurethane ¹ | | |
| | Uncuffed Polyurethane | Uncuffed Silicone |
| Exit Site Infections | 45% | 10% |
| Thrombosis | 20% | 2.6% |

Our silicone Hemo-Cath® devices provide a biocompatible and flexible design, which may reduce kinking and vessel damage.

Available in a full range of device configurations, the Hemo-Cath provides vascular access for long term or short term hemodialysis for adult patients.

Features & Benefits

- **Material**
Catheters designed with biocompatible and flexible silicone material.
- **Longer Lengths**
Enable adequate treatment with tip placement in the IVC via the femoral approach.
- **Hemo-Cath® Hub Clip**
Curves the extensions for Subclavian or IJ insertions.
- **Fast Track Stylet**
Provides adequate rigidity during insertion.
- **MRI Safe**

