

ACUTE
DIALYSIS
CARE

trio-ct™

TRIPLE LUMEN CATHETER

TRIO-CT™ CATHETER TIP



*Three lumens...
One catheter...
Perfect Harmony*

TekMed™
1300 720 727
www.tekmed.com.au



FS521104



TekMed NZ™
0800 456 217
www.tekmed.co.nz

medCOMP®



Third lumen provides additional access for multiple therapies.

ACUTE
DIALYSIS
CARE


trio-ct™

TRIPLE LUMEN CATHETER

Power Injectable Lumen

Reduces the need for additional access line

Reduces the need for placement of an additional central catheter for fluid and medication administration



International


ORDERING INFORMATION

CATHETER SIZE	INSERTION LENGTH	CATHETER TYPE	BASIC SET 5/BOX
13.5F	12cm	Straight	AS01135121E.
13.5F	15cm	Straight	AS01135151E.
13.5F	20cm	Straight	AS01135201E.
13.5F	24cm	Straight	AS01135241E.
13.5F	30cm	Straight	AS01135301E.

Catheter Only Contents: (1) Catheter (1) Dilator (1) Attachable Suture Wing (3) End Caps (1) Patient ID Card (1) Chart Sticker

Basic Set Contents: (1) Catheter (1) Introducer Needle (1) Guidewire (1) Dilator (1) Attachable Suture Wing (3) End Caps (1) Patient ID Card (1) Chart Sticker

TRIO-CT™ BASIC SET (INTERNATIONAL)



EC REP

MPS Medical Product Service GmbH
Borngasse 20
35619 Braunfels
Germany

Medical Components, Inc.
1499 Delp Dr., Harleysville, PA 19438, U.S.A.
P: 215.256.4201 F: 215.256.1787
www.medcompnet.com

CE2797
PN2653 Rev.A 5/19