



TRI-FLOW

Indications for Use: The Medcomp® Tri-Flow Catheter is indicated for use in attaining Short-Term vascular access for Hemodialysis and Apheresis. It may be inserted percutaneously and is primarily placed in the internal jugular vein of an adult patient. Alternate insertion sites include subclavian vein or femoral vein as required. The Medcomp® Tri-Flow Catheter is intended to be used less than (30) days.

Contraindications: This catheter is intended for Short-Term vascular access only and should not be used for any purpose other than indicated in these instructions. Do not insert catheter in thrombosed vessels.

Refer to Instructions for Use provided with the product for complete instructions, warnings, precautions, and contraindications. Observe all instructions for use prior to using products. Failure to do so may result in patient complications.



Polyurethane Material	Tapered Tip Tip Design	11.5F French Size	Straight, Curved Extensions Configurations
---------------------------------	----------------------------------	-----------------------------	--



TRIPLE INTERNAL
LUMEN DESIGN



BASIC SET		
XTP3114MTE	11.5F X 12CM STRAIGHT TRI-FLOW CATHETER SET	10/BOX
XTP3116MTE	11.5F X 15CM STRAIGHT TRI-FLOW CATHETER SET	10/BOX
XTP3118MTE	11.5F X 20CM STRAIGHT TRI-FLOW CATHETER SET	10/BOX
XTP3119MTE	11.5F X 24CM STRAIGHT TRI-FLOW CATHETER SET	10/BOX

BASIC SET CONTAINS:

- (1) Catheter
- (1) Dilator
- (1) Guidewire
- (1) Introducer Needle
- (3) End Caps

BASIC SET - CURVED EXTENSIONS		
XTP3116IUSE	11.5F X 15CM CE TRI-FLOW CATHETER SET	10/BOX
XTP3118IUSE	11.5F X 20CM CE TRI-FLOW CATHETER SET	10/BOX

TekMed
1300 720 727
www.tekmed.com.au

TekMed NZ
0800 456 217
www.tekmed.co.nz



FS521104



Refer to the Table of Contents for Important Risk Information regarding this device.