



BasicLine2

Technology for people

Form and function.

For the best treatment requirements.





BasicLine2

Technology for people

BasicLine2 was specially developed to meet the needs of dialysis centers, oncology day units, infusion and transfusion medicine.

BasicLine2 improves comfort and well-being for patients and enables nursing staff to work ergonomically.

QUALITY
MADE IN
GERMANY



With the BasicLine2, **2 motors** ensure quick and smooth adjustment of the **back, seat and leg sections**. The segments are adjusted simultaneously.



The illuminated handset impresses with its instinctive operation.



The large foot switch enables quick and easy adjustment to the head-down position.



Electronics are designed for minimal maintenance and ease of servicing.

The comfortable head cushion ensures both a pleasant lying and sitting position.



The adjustable footrest helps to minimize discomfort.

Equipotential bonding increases patient safety.

Give your surroundings more colour

The **high-quality design** of the therapy table with your choice of upholstery color makes it an ideal addition to any hospital, clinic or doctor's surgery.

A **wide range of faux leather colors** also helps you to create a pleasant and relaxing atmosphere in your environment.

The **pleasant feel and excellent material properties** ensure long-lasting, easy-care and comfortable use.

Hygiene in perfection

Smooth surfaces and a three-part lying surface ensure quick and easy cleaning - even in the spaces in between. The absence of Velcro as a fastening method also significantly enhances hygiene conditions.

The **upholstery** is also permanently protected against staining with a **PVC protective cover**. The protective cover is available for the leg section, the armrest and the headrest cushion.



Innovative overload protection

Thanks to the **innovative overload protection system**, the armrest springs back when the load exceeds 25 kg, thus **increasing safety for patients** when getting into the chair. At the same time, the overload protection prevents damage reducing repair costs and downtime. **It can be moved to the desired position with just one lever position and rotate completely** so that staff can use it instinctively.



Technical data

Dimensions

Width incl. armrests	90 cm
Width of the seat	58 cm
Overall length	203 – 230 cm
Armrests	52 x 15,5 cm
Total weight	68 kg

Adjustment options

Backrest angle	75° – 20°
Seat section angle	30° – 5°
Leg section angle	5° – 30°
Footrest	30 cm

Safety

Safe working load	175 kg
-------------------	--------



Manufacturer:
MEDICAL GmbH Medizintechnik
Ravensberger Straße 65
32278 Kirchlingern, Germany



Distributed by:
MEDICAL Sales GmbH & Co. KG
Wilhelmstraße 14 A
61381 Friedrichsdorf, Germany

Individual accessories



Paper roll holder



Handle on backrest



IV-pole



Cover leg segment



Console table



Adjustable footrest



USB port



Rechargeable battery



2nd hand switch

Note: The product shown in this brochure contains optional accessories that are not included in the standard version.



FS521104

