



Polyurethane
Material

Tapered Tip
Tip Design

14F
French Size

**Straight,
Pre-Curved**
Configurations



COAXIAL INTERNAL LUMEN DESIGN
Round venous lumen provides greater hemodynamic capabilities.



TAPERED TIP
Reduces vessel trauma and provides an easier insertion.

360° ARTERIAL PORTS
Help eliminate vessel occlusion and arterial insufficiency.



	BASIC SET	
DFXL144MTE	14F X 12CM 400XL, STRAIGHT	5/BOX
DFXL146MTE	14F X 15CM 400XL, STRAIGHT	5/BOX
DFXL148MTE	14F X 20CM 400XL, STRAIGHT	5/BOX
DFXL149MTE	14F X 24CM 400XL, STRAIGHT	5/BOX
DFXL144IJSE	14F X 12CM 400XL, PRE-CURVED	5/BOX
DFXL146IJSE	14F X 15CM 400XL, PRE-CURVED	5/BOX
DFXL148IJSE	14F X 20CM 400XL, PRE-CURVED	5/BOX
DFXL149IJSE	14F X 24CM 400XL, PRE-CURVED	5/BOX

BASIC SET
CONTAINS:
(1) Catheter
(1) Dilator
(1) Guidewire
(1) Introducer Needle
(2) Injection Caps
US and Foreign Patents Pending.



TekMed™

1300 720 727

12 Bellevue Crescent Preston Victoria 3072 Australia

T (03) 9487 9999 F (03) 9484 8792

info@tekmed.com.au www.tekmed.com.au



F5521104