




DOUBLE "D" INTERNAL LUMEN DESIGN

Polyurethane
Material

Tapered Tip
Tip Design

9F, 11F, 12F
French Sizes

Straight,
Curved
Extensions
Configurations



TAPERED TIP

Reduces vessel trauma and provides an easier insertion.

THERMOSENSITIVE MATERIAL

Softens to the body's temperature.

FIVE LARGE ARTERIAL HOLES

Strategically placed to prevent loss of blood flow.

BASIC SET		
XTP94MT=	9F X 12CM DOUBLE LUMEN CATHETER SET	10/BOX
XTP96MT=	9F X 15CM DOUBLE LUMEN CATHETER SET	10/BOX
XTP98MT=	9F X 20CM DOUBLE LUMEN CATHETER SET	10/BOX
XTP94IJS=	9F X 12CM IJ DOUBLE LUMEN CATHETER SET	10/BOX
XTP96IJS=	9F X 15CM IJ DOUBLE LUMEN CATHETER SET	10/BOX
XTP114MT=	11F X 12CM DOUBLE LUMEN CATHETER SET	10/BOX
XTP116MT=	11F X 15CM DOUBLE LUMEN CATHETER SET	10/BOX
XTP118MT=	11F X 20CM DOUBLE LUMEN CATHETER SET	10/BOX
XTP119MT=	11F X 24CM DOUBLE LUMEN CATHETER SET	10/BOX
XTP114IJS=	11F X 12CM IJ DOUBLE LUMEN CATHETER SET	10/BOX
XTP115IJS=	11F X 13.5CM IJ DOUBLE LUMEN CATHETER SET	10/BOX
XTP116IJS=	11F X 15CM IJ DOUBLE LUMEN CATHETER SET	10/BOX
XTP118IJS=	11F X 20CM IJ DOUBLE LUMEN CATHETER SET	10/BOX
XTP126MT=	12F X 15CM DOUBLE LUMEN CATHETER SET	10/BOX
XTP128MT=	12F X 20CM DOUBLE LUMEN CATHETER SET	10/BOX
XTP129MT=	12F X 24CM DOUBLE LUMEN CATHETER SET	10/BOX
XTP125IJS=	12F X 13CM IJ DOUBLE LUMEN CATHETER SET	10/BOX
XTP126IJS=	12F X 15CM IJ DOUBLE LUMEN CATHETER SET	10/BOX
XTP128IJS=	12F X 20CM IJ DOUBLE LUMEN CATHETER SET	10/BOX

BASIC SET CONTAINS:

- (1) Catheter
- (1) Dilator
- (1) Guidewire
- (1) Introducer Needle
- (2) Injection Caps

CATHETER ONLY CONTAINS:

- (1) Catheter
- (2) Injection Caps

CATHETER ONLY		
XTP114CT=	11F X 12CM DOUBLE LUMEN CATHETER	10/BOX
XTP116CT=	11F X 15CM DOUBLE LUMEN CATHETER	10/BOX
XTP118CT=	11F X 20CM DOUBLE LUMEN CATHETER	10/BOX
XTP114IJC=	11F X 12CM IJ DOUBLE LUMEN CATHETER	10/BOX
XTP116IJC=	11F X 15CM IJ DOUBLE LUMEN CATHETER	10/BOX
XTP118IJC=	11F X 20CM IJ DOUBLE LUMEN CATHETER	10/BOX



TekMed™

1300 720 727

12 Bellevue Crescent Preston Victoria 3072 Australia

T (03) 9487 9999 F (03) 9484 8792

info@tekmed.com.au www.tekmed.com.au

